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
***CARBON FREE INVESTMENT  
FOR PARTIAL DENTURES,  
PRESSABLE AND CERAMIC***

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### Carbon Free Investment for partial dentures, pressable and ceramic

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#### Tech support

Working time:	4 minutes
Settling expansion:	0.6%
Thermal expansion:	1.4%
Burnout Range:	1300° - 1800° F (704°C - 982°C)

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#### Burn-out and Casting

Conventional Place molds into furnace at room temperature and do a conventional overnight burnout. Soak mold at 1800°F (982°C) for 60 minutes and cast. No maximum number of rings.

Recommended Digital Approach - Semi-Rapid\* For 2 rings: Place molds into furnace at room temperature and ramp up to 1800°F (940°C) over an hour. Soak molds for 75 minutes and cast. Add 8 minutes for each additional ring up to 4. Make sure sprue hole is cherry red prior to casting.

Rapid: place molds into furnace at 1100°F (595°C). Raise temperature to 1725°F (982°C) over 30 minutes. Maximum 12 rings in furnace

Polvo / Powder	Líquido Total / Total Liquid
(gramos/grams)	(ml)
620	155
620	150
500	125
400	100
100	25

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#### WARNING -

This product contains free silica. Inhalation of dust may be harmful to your health and excessive inhalation will increase the risk of serious respiratory disease (silicosis). Take appropriate measures to avoid breathing dust. Immediate clean-up of spills is highly recommended. This product contains chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

#### WARNING -

Causes eye and skin irritation. Harmful or fatal if swallowed. May be absorbed through the skin. Inhalation may cause respiratory irritation. Prolonged overexposure may cause lung, liver and kidney damage. Contains ethylene glycol, which may cause birth defects based on animal studies.

**First Aid Procedure If inhaled:** Remove to fresh air. If breathing is difficult, have qualified personnel administer oxygen. Get immediate medical attention. If in eyes: Immediately flush eyes with plenty of water. Get immediate medical attention. If on skin: Immediately wash skin thoroughly with soap and water. If rash, irritation or symptoms develop: Get medical attention. In case of fire: Use carbon dioxide, alcohol foam or dry chemical.

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KEEP CONTAINER TIGHTLY CLOSED  
STORE IN A COOL, DRY PLACE  
FOR DENTAL LAB USE ONLY