

NEW

MODELO

DENTAL STONE PLASTER



1. WORK MODELS
2. STUDY MODELS
3. ANTAGONISTS

4. FLAKS
5. BASE

MODELO

DENTAL STONE PLASTER

Stone plaster resistance type III.

Manufactured under international standards according to the ANSI/ADA No.25 Dental Gypsum Products and ISO 6873. Guaranteeing a reliable, smooth, and stable duplicate for any required dental procedure.


Uses:

Work Models
Study Models
Antagonists
Flaks
Base

Advantages:



Smooth finishing surface



Compatibility with other filling materials



Dimensional stability



Hermetic bag to avoid unwanted hydration



Economic

Characteristics



Setting expansion
Resistance to compression
Regular Setting Time
Ratio water/powder

Specifications for Normal Setting

0.10%
26.5 MPa (270 kg/cm ²) *
8'00" – 10'00" *
24 ml / 100 g

*May vary depending on the temperature and humidity from the ambience.


Presentations:



White



Yellow



Blue

AENOR

GESTIÓN DE LA CALIDAD
ISO 9001

ER-0472/2016



Headquarter
Industria del Plástico 2113
Fracc. Industrial Zapopan, Norte
Zapopan, Jalisco, México.
mdc@dental.com.mx
+52 (33) 3833-0333

/mdcdentaloficial

Mexico, City.
Alzate 118, Eje 1 Norte
Col. Santa María de la Rivera
Delegación Cuauhtémoc
mdcmex@dental.com.mx
+52 (55) 5547-4747

www.dental.com.mx

USA
17800 S. Main St. 201
Gardena CA, USA 90248
customerservice@mdcdental.com
011 (310) 352 6980

01 (800) 3637-800